

Table of Contents

RPL™ (Reactive Pulsed Light) Technology-How does it work?	2
Intended Use	3
CONTRAINDICATIONS	4
WARNINGS	5
CAUTIONS	6
Possible Adverse Effects	7
Device description	10
User Instructions	12
Preliminary Test	12
Getting Started	13
Step 1: Prepare your skin	13
Step 2: Device operation	13
Cleaning & Storage	17
Troubleshooting	17
Customer Service	18
Device specification parameters	19
Labeling	20

RPL™ (Reactive Pulsed Light) Technology—How does it work?

RPL™ (Reactive Pulsed Light) is an advanced technology which uses clinically proven IPL (Intense Pulsed Light) technology in conjunction with a skin tone and skin contact sensor. The light emitted from the device penetrates the skin and is then absorbed by the melanin pigment in the hair shaft. The heat produced by this light goes into the follicle and generates a local thermal effect, also known as selective photothermolysis. The thermal effect disables the hair follicle which causes a reduction in hair growth while leaving the surrounding skin unaffected and cool.


Long-term hair reduction—Hair growth cycle

Two factors play an important role in optimal long-term hair reduction:

- **Hair growth cycle.** Hair growth involves three phases: Anagen, Catagen, and Telogen. Anagen is the growth phase, while

Catagen and Telogen are the resting phases. Permanent hair reduction can only occur in the Anagen phase, the phase of growth. Since hair growth is not synchronized, it is recommended to complete the full treatment protocol to ensure every hair will be exposed to treatment during its Anagen phase.

- **Skin tones (according to the Fitzpatrick scale).** Darker skin contains more melanin, which competes with the targeted hair for energy absorption. When more energy is absorbed into the skin, there is a greater risk of adverse effects.





Please make sure to carefully read this user guide before you start using the device. The user guide can also be found at: www.sensica.com

The device is intended for use on skin tone types I to IV according to the Fitzpatrick skin chart below (colored chart available in section 02 on the cover). Do not use on dark skin such as type VI. The use of IPL devices on naturally dark skin is forbidden.

Intended Use

Sensilight Mini device is an over-the-counter device intended for the removal of unwanted hair. The Sensilight Mini is also intended for permanent reduction in hair regrowth, defined as a long-term, stable reduction in the number of hair regrowing when measured at 6, 9, and 12 months after the completion of treatment regimen.

This device is intended for use on the entire body (e.g. legs, arms, underarms, bikini line, back, shoulders, chest and stomach) and on permitted facial areas (i.e. above upper lip, cheek, jaw, chin and neck).

Note: This device is especially designed for your own personal treatment and is intended for a single user and home use only.

**Skin
Tones**



This device works best on dark hair or hair that contains more melanin. Brown and light brown hair will respond, but typically will require more treatment sessions.



Never use around or near your eyes

Note: Hair removal by lasers or Intense Pulsed Light can cause increased hair growth in some individuals. Based upon current data, the highest risk groups for this response are females of Mediterranean, Middle Eastern, and South Asian heritage, treated on the face and neck.

Who should NOT use the device?

Carefully read the information below before you begin a treatment with the device. If you still have any doubts or are unsure that the device is safe for your use, please consult with your physician.

CONTRAINDICATIONS

Do NOT use the device if you:

- Have an active implant such as a pacemaker, insulin pump or similar devices, since it may cause interference.
- Are pregnant, have a chance of being pregnant or are nursing since the device was not tested on pregnant or nursing women and may stimulate pigmentation.
- Have skin cancer, areas at risk of malignancy, or if you have received chemotherapy/radiation therapy within the past 3 months since the safety of such usage has not been established and the immune system can be compromised, which can lead to serious injury.
- Suffer from epilepsy, as using the device could trigger a seizure.
- Have damaged or broken skin, or a history of Herpes or Psoriasis outbreak in the treatment area since such usage can result in serious side effects.
- Suffer from keloid scar formation (growth of extra scar tissue where the skin has healed after an injury) in the treatment area since such usage has not been established and may result in serious skin complications.
- Have been taking Tetracycline antibiotics or medication with the following active ingredients within the last two weeks: Retinoid acid, Retina-A, Vitamin A (such as Accutane), as they are considered to be photosensitizing drugs and may result in skin injury.
- Are photosensitive or have other diseases related to light sensitivity, as device usage under these conditions may lead to serious injury.
- Suffer from severe or decompensated metabolic

conditions, such as diabetes, since it may cause a greater risk of delayed healing and infection in case of an adverse effect.

- Have had a skin peel or other skin resurfacing procedures within the last 8 weeks in the treatment area as, under these circumstances, device usage may result in stimulation of pigmentation.

Note: If you have undergone any other aesthetic procedures or have any uncertainty regarding potential side effects, consult with your physician before using the device.

WARNINGS

Do NOT use the device:

- On naturally dark skin (colored chart available in section 02 on the cover).
- If you have dark brown or black spots such as: around the nipples, freckles, birth marks, moles, or blisters on the area to be treated.
- On mucous membranes such as: nostrils, lips of the mouth, eyelids, ears, genital

area and the anus, as it may lead to serious injury.

- Over piercings or any metal items such as earrings or jewelry on the area to be treated, as it may cause burns. Over areas with tattoos or permanent make-up, cover tattoos or permanent make-up with a white plaster before treatment. Do not use the device directly over the tattoo or permanent make-up, even when covered.
- **On tanned skin or after sun exposure as this can cause serious burns or skin injury. Avoid more than 15 minutes of sun exposure for 4 weeks before your treatment and for at least 2 weeks post treatment. If during this period your treated skin area is exposed to sunlight, make sure to apply SPF 30 or higher sunscreens.**

Additionally:

- Never look directly at the light coming from the applicator and lamp, as it may cause eye injury.
- Never use flammable liquids such as alcohol (including perfumes, sanitizers, or other applications containing alcohol) or acetone on the skin one hour before using the device as it may cause burning of the skin.
- Do not treat the same area of skin more than once per treatment, as it may cause adverse effects.

CAUTIONS

- KEEP OUT OF REACH OF CHILDREN.
- This device is not intended for use by individuals with reduced physical, sensory or mental capabilities.
- Do not place, store or use this device near bathtubs, showers, basins or other vessels containing water, as it may lead to severe electrocution.
- Do not open or repair the device, as it may expose you to dangerous electrical components

which may cause serious bodily damage.

- Do not leave the device unattended while plugged into the electricity, as it is a risk of electrocution.
- Do not use the device if any part of the device or its accessories are: not working properly, appear damaged, broken, loose, missing or wet (take extra care since UV and IR are emitted from this device), or if you see or smell smoke.
- Before using this device, make sure it is clean, dry and intact. Cleaning instructions are detailed under the **Cleaning & Storage** section.



Use the device only with the supplied external power and accessories. For your safety, never use DC adapters which are not approved by the manufacturer.

Possible Adverse Effects

When using this device according to the Instructions for Use, adverse effects are rare. However, every aesthetic procedure may cause some degree of adverse effects (as listed in the table below). In case adverse effects appear, stop using the device, and consult with your physician if necessary.

Mild adverse effects should be spontaneously relieved within a few hours.

Severity of Adverse Effect	Frequency of Adverse Effect	Adverse Effects
Mild	Common	Skin Irritation Skin Redness Skin Sensitivity
Moderate	Rare	Skin Burns Scarring Pigment Changes Excessive Redness and Swelling
Severe	Extremely Rare	Infection Bruising

Any serious incident that has occurred in relation to the device should be reported to us (via customer service).

For your convenience, additional explanations are given below.

- **Skin Irritation** - The use of Pulsed Light Hair Removal System is generally very well- tolerated, the majority of users do feel a sensation of warmth and a mild stinging sensation during treatment or a few minutes after. Anything beyond this skin irritation – discomfort is abnormal and means that either the energy level is too high, or you shouldn't continue using the device as you are unable to tolerate the hair reduction application.
- **Skin Redness** - Redness on the treated area may appear right after using the application or within 24 hours of using the device. Consult your physician in case the redness does not disappear within 2 to 3 days.
- **Skin Sensitivity** - You may feel that your skin in the treated area is more sensitive so you may encounter dryness or flaking of the skin.
- **Skin Burns** - Rarely, burns on the treated area can occur following the use of the hair removal system. The burns can require a few weeks to heal and, extremely rarely, may leave a visible permanent scar.
- **Scarring** - In rare cases, permanent scarring may occur. Usually when scarring occurs it is in the form of a flat and white lesion on the skin (hypotrophic). However, it can be large and red (hypertrophic) or

large and extend beyond the margins of the injury itself (keloid). Subsequent aesthetic treatments may be required to improve the appearance of the scar.

- **Pigment Changes** – Use of the device may cause temporary skin hyperpigmentation (increased pigment or brown discoloration) or hypopigmentation (whitening) to the surrounding skin. The risk is higher for people with darker skin tones. Usually discoloration or changes to skin pigment are temporary and permanent hyperpigmentation or hypopigmentation rarely occur.
- **Excessive Redness and Swelling** – Very rarely, the treated area skin may become very red and swollen. This side effect is very common in sensitive areas of the body. The redness and swelling should be decreased after 2 to 7 days and you should treat the area by applying ice. You should avoid exposure to sun.
- **Infection** – In rare cases, infection of the skin is a possible risk following a skin burn caused by the application.
- **Bruising** – Blue-purple bruising is exceedingly rare by using the device and can last 5 to 10 days. As the bruising fades, there may be a rust-brown discoloration of the skin (hyperpigmentation) that can be permanent.

Device description

- 1 ON/OFF button and energy level button
- 2 Trigger button
- 3 Status indicator:
 - Lack of contact
 - Non-allowed skin tone
 - Device not usable
- 4 1-3 energy level display
- 5 Skin tone and contact sensor (Fitzpatrick sensor)
- 6 Power supply socket
- 7 Treatment window
- 8 Lamp



User interface description

- 1. ON/OFF and energy selection button** - Pressing the button (after the device is connected to the power supply) will turn the device ON at energy level 1. Each press will change the energy, from the lower level to the higher, for a total of 3 different energy levels. Select the energy level according to what you found most comfortable in the preliminary test. The energy level display indicators will light up green according to the selected energy level (e.g. at energy level 1, one green light will be lit, at energy level 2, two green lights will be lit). The device will switch to standby mode when left unattended for 3 minutes or by pressing the ON/OFF button for 3 consecutive seconds.

2. Trigger button - (flash pulse activation) marked as "sun" on the upper side of the device. The trigger button will activate the flash lamp if the treatment window is in full contact with the skin, and the skin tone is in within the allowed range. Pressing the button once will activate a single flash, pressing and holding the button will activate flashes continuously as long as the button is pressed.

Status Indicator	Status	Description
No light	OFF	Device not working (not connected to the power supply).
Steady green light	ON	Flash is ready to be delivered.
Blinking green light	Standby	Device is in standby mode (before powering on or when left unattended for 3 minutes).
Momentary orange light	Flash blocked	The device is not in full contact with the skin – no flash will be delivered.

User Instructions

The device is suitable for use on: legs, arms, underarms bikini line, back, shoulders, chest, stomach, and permitted facial areas (above upper lip, cheek, jaw, chin and neck).

What to expect

During your treatment, it is normal to experience:

- A flash of light with each press on the trigger button – the emitted light will not harm your eyes.
- A pop noise with each flash of light.
- A mild, stinging sensation.

Preliminary Test

The preliminary test is intended to determine your skin's reaction to treatment with the device and help you set the optimal energy level for each treatment area.

Perform a preliminary test on a hidden patch of skin for the area you intend to treat 48 hours prior to the first full treatment.

Due to varying levels of sensitivity, body areas may require different energy levels; this will be determined by separate preliminary tests.

Note: Do not treat full areas with energy levels that you have not tested before!

Do not treat the tested area for at least 48 hours after the preliminary test has been carried out.

1. Prepare your skin as described below in "Step 1: Prepare your skin". Instructions for operating the device can be found below in "**Step 2: Device operation**".
2. Select the lowest energy level and deliver a flash.
3. You can increase the energy level as long as it feels comfortable. You may continue to do this up to the maximum energy level unless you begin to experience discomfort.
4. If you do not feel comfortable or see adverse effects such as severe redness or blistering stop the test immediately.
5. Wait 48 hours and examine the tested area:
6. If the skin in the test area seems normal (there is no reaction or just a mild

redness), you can proceed with a full treatment on the highest comfortable energy level that you tested.

Getting Started

STEP 1: PREPARE YOUR SKIN

1. Remove all jewelry.
2. Make sure the hair on your skin in the treatment area is shaved or trimmed (between treatments you can shave, but don't wax or pluck hair), and that the skin is free of any residue from: powders, creams, makeup, antiperspirants or deodorants.
3. Never use flammable liquids such as alcohol (including perfumes, sanitizers, or other applications containing alcohol) or acetone on the skin for one hour before using the device, as it may cause burning of the skin.

Important: Do not use any cream or gel while using the device since it may prevent normal operation.

Remember: The device is equipped with a built-in skin-tone sensor, allowing the device to work only on suitable skin tones. This unique safety feature will prevent a flash from occurring on specific areas where your skin is too dark or too tanned, even if your general skin type is lighter.

STEP 2: DEVICE OPERATION

Note: A preliminary test should be conducted before proceeding with a full treatment. See "Preliminary Test" section above. Make sure the device is clean, dry and intact.

1. Turn the device on. When the device is turned on it will be automatically set to the lowest energy level and the energy level display will show one green light. Manually change energy levels by pressing the "ON/OFF" and energy selection button. You may decrease or increase the energy level according to your comfort level, as determined by the preliminary test.

- Place the treatment window on the desired treatment area, ensuring full contact with the skin. As a safety precaution the device is equipped with a safety sensor that will not trigger a flash if:
 - The device is not properly applied to the skin. In this case, the status indicator will turn "momentarily" orange.
 - Your skin tone is too dark/tanned. In this case, the status indicator will turn "momentarily" red.
- Apply the treatment window firmly against your skin and press the trigger button to release a flash. Following each flash, move the device to the next spot. After each flash, the device will take 1-3 seconds to be ready for the next one.
- For continuous flashes, keep the trigger button pressed and start moving the device with each flash to the next spot without detaching it from the skin. The device will immediately recharge and will be ready for the next flash. Repeat process until the entire treatment area has been covered.
- After the treatment, turn off the device by pressing the "ON/OFF" button. Unplug the unit after each treatment session.
- Perform treatment sessions according to the "**Recommended Treatment**" section.

Important note: Try not to overlap pulses – repeated pulses on the same spot may cause burns.

What to expect:

Different individuals and different treatment areas may not react equally to RPL technology.

You should begin to notice a reduction in hair growth by the 6th treatment. Once the treatment protocol of 12 treatments is complete, you should see further hair growth reduction. Any hairs that remain should be finer and lighter in color. In order to maintain results, we recommend using the device once a month per treatment area.

Important: Facial treatments should be performed at the medium or low energy level.

Notes:

- It is recommended to use the highest energy setting that you found comfortable during the preliminary test.
- Lighter hair colors may require more treatments in order to see results.
- The hair will continue to grow for up to two weeks after a hair reduction session. However, as you progress through the full treatment protocol, hair will naturally and gradually fall out. The time until initial results can be noticed varies between users and depends on factors such as: age, hair color and skin tone.
- The number of sessions required for achieving long term results differs for every treatment area and between different users.
- This device is suitable for men. However, male body hair, especially on the chest, will require more treatment sessions.
- Maintenance will be needed from time to time if growth is still visible.

Recommended treatment protocol:

Treatment areas depending on hair color and their treatment protocol:

Hair Color	Fitzpatrick Skin Level	Energy Levels	Treatment Area	Average Number of Treatments	Treatment Interval
Light Brown to Brown	I–IV (Inclusive)	1–3	<ul style="list-style-type: none">• Legs• Arms• Underarms• Bikini Lines	10–12	The 1–4 treatments will be approximately two weeks apart, followed by four-week time intervals.
Brown to Black	I–IV (Inclusive)	1–3	<ul style="list-style-type: none">• Legs• Arms• Underarms• Bikini Lines	8–10	

Important: Facial treatment should be performed at 1–4 energy levels only.

Post-treatment care

- After treatment, you may use a calming lotion or try Sensica's post RPL lotion (available at www.sensica.com).
- Avoid using peeling and bleaching creams or similar products for at least 24 hours, as this irritates the treated skin.
- Avoid direct sun exposure for at least 2 weeks post treatment. If exposing recently treated skin to sunlight during this time period, apply sunscreen of SPF 30 or higher.

Cleaning & Storage

1. Make sure to turn off and unplug the device before cleaning and storage.
2. It is recommended that the treatment window area be cleaned after every session by using a dry and clean cloth.
3. If there is dirt or debris on the treatment window, clean the treatment window. If you can't get the device clean, contact our Customer Service.

Troubleshooting

If the device does not turn on:

- Make sure that the power supply is plugged properly into the device and the electrical outlet.
- Make sure the power outlet has power.

If the device does not trigger a pulse:

- Verify there is full contact between the device and your skin. The trigger button will activate a pulse only when the device is in full contact with the skin.
- Verify that the area you are trying to treat is not too small, curved or bony. The entire treatment window should touch the treatment area, pulling the skin flat can help to achieve this.
- If the skin tone identified by the skin tone and contact sensor is too dark, the flash activation will be blocked, and the status indicator will turn "momentarily" red.

In case of general error:

- Reset the device by unplugging the power cord from the electrical outlet on the wall, wait 10 seconds


and then reconnect it. Then try turning on the device after resetting.


- In case of a pulse delay after repositioning the device on a different treatment area, you should consider performing the treatment in a cooler environment. The device may pulse at a slower rate in a warm environment.
- In case of dirt or debris is stuck on the treatment window, you should clean the treatment window according to the cleaning section. If you can't get the area clean, contact our customer service.

If the problem is not solved, please contact our customer service

Customer Service

For more information about the device, if your device is broken, stopped working, or if you need more assistance, please contact us:

 855-SENSICA
(855-736-7422)

 support@sensica.com

Device specification parameters

Manufactured by:

ELG trade Ltd.

Emitted light spectrum:

475 [nm] to 1,200 [nm] max.

Emitted energy flux:

5 [J/cm²]

Selectable energy levels:

3 levels

Treatment area (window):

3 [cm²]

Speed:

1 flash every 3 seconds

External Power Supply:

12[V], 1.8[A]

100–240 [VAC], 50–60[Hz]

Device weight and size:

220 [gr], 135 X 45 X 75 [mm]

7.8 [Oz], 5.3 X 1.77 X 2.96 [in]

Technology:

Home use IPL

Storage Instructions

Storage and transport conditions

Temperature	14°F to 158°F (-10°C to 70°C)
--------------------	----------------------------------

Humidity	131°F (0 to 90% RH at 55°C)
-----------------	--------------------------------

Pressure	non-condensing
-----------------	----------------

Operational conditions

Temperature	50°F to 104°F (10°C to 40°C)
--------------------	---------------------------------

Humidity	30% to 70% RH
-----------------	---------------

Max operating altitude	10,000 ft.
-------------------------------	------------

Labeling

Model: IPL030US



Double insulation



Type BF equipment



Consult Instructions for use



Waste Electronic and Electrical Equipment



Keep away from eyes



Please read all warnings and safety precautions prior to use



Do not use the device on tanned skin or after sun exposure



Do not let the device become wet. Keep away from bathtubs, showers, basins, or other vessels containing water: Keep dry



Use the device only with the supplied external power and accessories. For your safety, never use DC adapters which are not approved by the manufacturer.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

1. This device may not cause harmful interference.
2. This device must accept any interference received, including interference that may cause undesired operation.

Electromagnetic compatibility (EMC) and Specifications

- The Sensilight Mini needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided.
- Portable and mobile RF communications equipment can affect Sensilight Mini.
- The use of accessories or replacement parts, other than those sold/specified by sensica International, may result in increased emissions or decreased immunity of the Sensilight Mini.
- The Sensilight Mini should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, the Sensilight Mini should be observed to verify normal operation in the configuration in which it will be used.

EMC Information

Guidance and manufacturer's declaration – electromagnetic emissions
as per IEC 60601-1-2: 2007 - Table 1

The Sensilight Mini is intended for use in the electromagnetic environment specified below. The customer or the user of the Sensilight Mini should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The Sensilight Mini uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Sensilight Mini is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration – electromagnetic immunity as per IEC 60601-1-2: 2007 - Table 2

The Sensilight Mini is intended for use in the electromagnetic environment specified below. The customer or the user of the Sensilight Mini should assure that it is used in such an environment.


Immunity Test	IEC 60601 Test Level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines	± 2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) *Not Applicable *± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_T (>95 % dip in U_T) for 0,5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 s	<5 % U_T (>95 % dip in U_T) for 0,5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Sensilight Mini requires continued operation during power mains interruptions, it is recommended that the Sensilight Mini be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE U_T is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration – electromagnetic immunity

As per IEC 60601-1-2: 2007 - Table 4

The Sensilight Mini is intended for use in the electromagnetic environment specified below. The customer or the user of the Sensilight Mini should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance level	Electromagnetic environment – guidance
<p>Conducted RF IEC61000-4-6 Radiated RF IEC 61000-4-3</p>	<p>3Vrms 150 kHz to 80 MHz 3 V/m 80 MHz to 2,5 GHz</p>	<p>3 V 3 V/m 80 MHz to 2,5 GHz</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the Sensilight Mini, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> <p>$d=3.5/3 \sqrt{P}=1.17\sqrt{P}$</p> <p>$d=3.5/3 \sqrt{P}=1.17\sqrt{P}$ 80 MHz to 800 MHz</p> <p>$d=7/3 \sqrt{P}=2.33\sqrt{P}$ 800 MHz to 2.5 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey a, should be less than the compliance level in each frequency range b.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol: </p>

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Sensilight Mini is used exceeds the applicable RF compliance level above, the Sensilight Mini should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Sensilight Mini.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the Sensilight Mini as per IEC 60601-1-2: 2007 - Table 6

The Sensilight Mini is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Sensilight Mini can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Sensilight Mini as recommended below, according to the maximum output power of the communications equipment.

Rated Maximum Output Power of Transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d=3.5/\sqrt{P}$ $\sqrt{P}=1.17\sqrt{P}$	80 MHz to 800 MHz $d=3.5/\sqrt{P}$ $\sqrt{P}=1.17\sqrt{P}$	800 MHz to 2,5 GHz $d=7/\sqrt{P}$ $\sqrt{P}=2.33\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.70	3.70	7.37
100	11.67	11.67	23.33





For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.




Treatment log

Treatment date	Treated area	Comments
----------------	--------------	----------



Treatment log

Treatment date	Treated area	Comments
----------------	--------------	----------



Treatment log

Treatment date	Treated area	Comments
----------------	--------------	----------

